




A Fully Guided Surgical Solution for the LOCATOR Implant System



The comprehensive Guided Surgical Kit from Zest Dental Solutions is designed to be used with the trusted LOCATOR® Implant System.

Our Guided Surgical Kit is built for performance, efficiency, precision, and peace of mind, giving clinicians the tools they need to drill and place implants confidently while exceeding patient expectations.

This all-in-one kit does not rely on keys or spoons, thereby maximizing efficiency, and includes the necessary instruments to support multiple surgical access options. A single guide sleeve is used for straightforward surgical guide creation and workflow. Three total drill lengths per size improve efficiency and treatment flexibility, and are supported by a color coded drill protocol. All of this in a compact, well organized setup helps to instill confidence and deliver an exceptional surgical experience.

Available from any participating Zest selling partner or directly from Zest Dental Solutions.

For more information go to ZestDent.com or call 800-262-2310.



FEATURES & BENEFITS

- **Simplicity by design.** Keyless guided surgical system utilizes a single guided sleeve to simplify instruments and surgical protocols.
- **Easy-to-follow, color-coded drill protocol.** Color scheme is clearly depicted on the kit's front for quick reference during surgery.
- **Small, compact footprint with angled presentation.** Fits in seamlessly with your operatory and allows for easy storage and access.
- **Stay Organized!** Store LOCATOR CORE tool and torque wrench in one location.
- **Peace of mind.** Drill depth guides include three drill lengths with lasermark line at 4mm short of prolongation value.
- **Built in accuracy.** Measurement tool provides a quick reference check for implants or drill lengths.
- **Precision-based flexibility.** Allows for multiple surgical approaches support (flapped & flapless) including NEW Starter drill/tissue punch combo instrument.



LEARN MORE AT
WWW.ZESTDENT.COM

Available from any participating Zest selling partner or directly from Zest Dental Solutions.

For more information go to ZestDent.com or call 800-262-2310.