



# IDEAL FOR BLOCKING OUT UNDERCUTS WHILE PROCESSING OVERDENTURES CHAIRSIDE






## BLOCK OR AVOID UNDERCUTS USING CHAIRSIDE® BLOCK OUT MATE-

Zest Dental Solutions continually listens to clinicians about how the company can design, manufacture and source products that make implant-retained overdenture treatment simpler and more efficient, as well as creating a comfortable and satisfying procedure for patients. The latest advancement in the CHAIRSIDE Product Portfolio supports that mission

- **Introducing CHAIRSIDE Block Out Material** from Zest Dental Solutions, **a light-cured protective material that removes undercuts while processing overdentures chairside.**




### PATIENT COMFORT

CHAIRSIDE Block Out Material is an acrylic resin that remains cool and comfortable during light curing, preventing discomfort and harm to the tissue.

### EASY HANDLING

Viscous material stays put using syringe delivery. Contrasting blue color is easy to see, flexible and simple to remove.



### MULTIPLE USES

**CHAIRSIDE Block Out is a superior material for blocking out undercuts for LOCATOR® and SATURNO® cases - it is also ideal to:**

- Protect tissue while bleaching, sandblasting or using air abrasion
- Seal around rubber dams
- Repair torn rubber dams
- Maintain reservoir space for bleaching trays

ITEM NUMBER	DESCRIPTION
09574	<ul style="list-style-type: none"><li>• Four (4) 1.2mL Syringes</li><li>• Twenty (20) Mixing Tips</li></ul>



**ZEST DENTAL SOLUTIONS**  
2875 LOKER AVENUE EAST  
CARLSBAD, CA 92010  
USA

TEL: 800.262.2310  
FAX: 760.743.7975  
EMAIL: [SALES@ZESTDENT.COM](mailto:SALES@ZESTDENT.COM)  
[WWW.ZESTDENT.COM](http://WWW.ZESTDENT.COM)