



LOCATOR ATTACHMENT SYSTEM INSERT RETENTION STRENGTH



	EXTRA LIGHT	LIGHT	REGULAR
LOCATOR® Standard Range Insert Up to 10° Angle Correction / 20° Between Two Implants			
LOCATOR® Extended Range Insert Up to 20° Angle Correction / 40° Between Two Implants			

LOCATOR R-TX ATTACHMENT SYSTEM INSERT RETENTION STRENGTH



	LOW	MEDIUM	HIGH
LOCATOR R-Tx® Standard Range Insert Up to 30° Angle Correction / 60° Between Two Implants			
LOCATOR R-Tx® Limited Range Insert Up to 5° Angle Correction / 10° Between Two Implants			