

LOCATOR R-TX[®] ABUTMENTS WITH LOCATOR[®] OVERDENTURE IMPLANTS NOW AVAILABLE!

NEW LOCATOR IMPLANT SYSTEM is now compatible with **LOCATOR R-TX ABUTMENTS** offering 25% greater wear resistance, an impressive 60° of total implant divergence and dual retentive surfaces for easier seating of the overdenture. LOCATOR R-Tx delivers greater durability with easier maintenance and is the perfect companion to the AFFORDABLE LOCATOR IMPLANT¹.



NEW






DURATEC®
 ABUTMENT
 COATING

30% INCREASE IN STRENGTH¹


25% GREATER WEAR RESISTANCE¹

**65% REDUCTION IN SURFACE
 ROUGHNESS¹**

Patented pivot technology that allows a total of **60° OF IMPLANT DIVERGENCE.**



Limited Range inserts are intended for patients with **10° OR LESS OF IMPLANT DIVERGENCE**, providing a tighter, more secure fit for patients wanting more stability.



ORDERING INFORMATION



2.4mm		2.9mm		3.5mm	
Part	Length	Part	Length	Part	Length
07381	10mm	07386	10mm	07501	8mm
07382	12mm	07387	12mm	07502	10mm
07383	14mm	07388	14mm	07503	12mm
				07504	14mm



3.9mm		4.4mm		4.9mm	
Part	Length	Part	Length	Part	Length
07505	8mm	07509	8mm	07513	8mm
07506	10mm	07510	10mm	07514	10mm
07507	12mm	07511	12mm	07515	12mm
07508	14mm	07512	14mm	07516	14mm



R-Tx Abutments for LOCATOR Implants	
Part	Description
32802-02-SB	2.5mm cuff height
32802-03-SB	3mm cuff height
32802-04-SB	4mm cuff height
32802-05-SB	5mm cuff height
32802-06-SB	6mm cuff height


IMPROVED PERFORMANCE

"I love having the option of being able to use LOCATOR R-Tx on LOCATOR Overdenture Implants. Since LOCATOR R-Tx abutments accommodate up to 60° it allows me to use LOCATOR implants in more challenging ridge situations. I've also noticed the new DuraTec coating on the R-Tx abutments helps to reduce plaque accumulation and helps to reduce maintenance over time."

-Michael Scherer, DMD, MS, FACP



LOCATOR R-Tx[®]



VISIT
ZESTDENT.COM
24 HOURS A DAY | 7 DAYS A WEEK

TO PLACE ORDERS ONLINE OR TO FIND THE
LATEST SUPPORT TOOLS FOR R-Tx[®] ABUTMENTS
WITH LOCATOR OVERDENTURE IMPLANTS

**AVAILABLE DIRECTLY FROM ZESTDENT.COM OR FROM
ONE OF OUR PARTICIPATING SELLING PARTNERS**



2875 LOKER AVENUE EAST
CARLSBAD, CA 92010
USA

TEL: 800.262.2310 | FAX: 760.743.7975
EMAIL: MARKETING@ZESTDENT.COM
WWW.ZESTDENT.COM

¹ Data on file, Zest Dental Solutions

©2022 Zest Anchors LLC. All rights reserved. DuraTec and Zest Dental Solutions are registered trademarks of ZEST IP Holdings, LLC.

L9165 REV B 08/2022