

INDICATIONS

The LOCATOR Fixed™ Attachment System is designed to support fixed, partial or full arch restorations on endosseous dental implants in the mandible or maxilla for the purpose of restoring masticatory function. It is used in fixed hybrid restorations that can be attached with a snap-in system.

IMPLANT COMPATIBILITY

Implant Manufacturer	Implant System	Implant Diameter
Ace Surgical	Infinity OCTAGON	3.3, 4.1, 4.8
	Infinity TRI-CAM	3.5, 4.3, 5.0
	Infinity Internal Hex	3.7, 4.1, 4.7, 5.1
Biohorizons	Tapered Internal	3.0
	Tapered Plus	3.8
	Laser-Lok	3.0
Biomet 3i	3iT3, NanoTite, OSSEOTITE	3.25, 3.75, 4.0
Blue Sky Bio	Quattro	3.3, 4.1, 4.8
	One Stage	3.3, 4.1, 4.8
Camlog	SCREW-LINE ROOT-LINE 2	3.3, 3.8, 4.3, 5.0
	SCREW-LINE	3.3, 3.8, 4.3, 5.0
Dentium	Superline	3.4, 3.8, 4.3, 4.8, 5.8
	Superline II	3.4, 3.8, 4.3, 4.8, 5.8
	Implantium	3.6, 3.8, 4.3, 4.8
Dentsply	Astra Tech OsseoSpeed EV	3.6, 4.2, 4.8
	Xive Frialit-2	3.4, 3.8, 4.5, 5.5
	Ankylos C/X	3.5
	Astra Tech OsseoSpeed TX	3.5, 4.0, 4.5, 5.0
DIO	UF(II)	3.8, 4.0, 4.5, 5.0, 5.5
Hiossen, Inc.	SS	3.5, 4.0, 4.5
	ET, TS	3.5, 4.0, 4.5, 5.0, 6.0, 7.0
	US	3.5, 3.3, 4.5, 4.0, 5.0, 6.0, 7.0
Implant Direct	Swish Plus	4.1, 4.8

Implant Manufacturer	Implant System	Implant Diameter
	SwishTapered	4.1, 4.8
	InterActive	3.2, 3.7, 4.3, 5.0
	ReActive	3.7, 4.2, 4.7, 5.7
	RePlant	3.5, 4.3, 5.0
	Legacy 1, 2, 3, 4	3.7, 4.2
	Legacy 2, 3, 4	3.2, 4.7, 5.2
	Legacy 1, 3	5.7
	Legacy 2, 4	5.7, 7.0
	Legacy 1	4.7
IDS	MegaGen Any Ridge	3.5, 4.0, 4.5, 5.0, 5.5, 6.0, 6.5, 7.0, 7.5, 8.0
IQ	Master	3.75, 4.2, 5.0, 6.0
	Smart	3.75, 4.2, 5.0, 6.0
KeyStone	PrimaConnex	3.5, 4.1, 5.0
	Genesis	3.5
	TILOBEMAXX	7.0, 8.0, 9.0
MIS	C1	3.3, 3.75, 4.2
	V3	3.3, 3.9, 4.3, 5.0
	SEVEN	3.3
	M4	3.3
Neodent	Helix GM	3.5, 3.75, 4.0, 5.0, 6.0
	Drive GM	3.5, 4.3, 5.0
	Titamax GM	3.5, 3.75, 4.0, 5.0
Nobel	Replace, Replace Select, NobelSpeedy	3.5, 4.3, 5.0
	NobelActive, NobelParallel CC, NobelReplace CC	3.5, 4.3, 5.5
	Branemark, Nobel Speedy, Groovy	3.3, 3.75, 4.0, 5.0
OCO Biomedical	Engage	3.25, 4.0, 5.0
Paltop	Advanced	3.25, 3.75, 4.2, 5.0
	Dynamic Conical	3.25, 3.75, 4.2, 5.0
	Conical Active	3.25, 3.75, 4.2, 5.0

Implant Manufacturer	Implant System	Implant Diameter
Southern Implants	Tri-Nex	3.5, 4.3, 5.0
	Tri-MAX7	7.0
Straumann	Roxolid SLActive, Roxolid SLA	3.3, 4.1, 4.8
	BLX	3.8
	TLX	3.5, 4.8, 6.5
Zest Dental Solutions	LOCATOR Overdenture Implant System	2.4, 2.9, 3.5, 3.9, 4.4, 4.9
Zimmer Dental	Trabecular Metal, Tapered Screw-Vent, Screw-Vent, Advent	3.7, 4.1, 4.7, 6.0
	Spline Reliance Cylinder	3.25, 4.0, 5.0
	Spline Twist	3.75